

BECCHISLINE
b e a u t y c o m p a n y

BL

Via Sommariva 31/1 - 10022 Carmagnola (TO)
Tel. +39 011.97.10.730 - info@becchisline.com
www.becchisline.com

Slimming, Lymphatic drainage

2 in 1 PRESSOTHERAPY

Snell System

Slimming method
without operator

MADE IN ITALY

www.becchisline.com

Body size and weight

Height: 10 cm

Width: 32 cm

Length: 41 cm

Weight: 5 kg

Size and overalls set

Height: 12 cm

Width: 46cm

Length: 91 cm

Weight: 6,6 kg

Technical features

Power supply: 230V / 50Hz

Electricity consumption: 350 W

Max pressure: adjustable from 30
to 120 mm hg

Electrical absorption: 160 VA

Security class: II / BF

Sternum protection: 2 fuses 1,6 A

Inflating sectors: 18 independent

Compression: 1 to 15 seconds per sector

General time: 5 to 60 minutes

N° of programs: 4

Infrared temperature: recommended from 45° to 50°

YOU CAN USE IT AS

IMPROVEMENT OF CIRCULATION

deflation promotes outflow
of extracellular impurities towards
the vessels by improving the
microcirculation lymphatic

LYMPHATIC DRAINAGE

has a remarkable draining function,
improves both blood circulation and
the lymphatic one

SLIMMING

Localized loss of cm

DRAINING

Remarkable draining function

