

**BECCHISLINE**  
b e a u t y c o m p a n y

**BL**

Via Sommariva 31/1 - 10022 Carmagnola (TO)  
Tel. +39 011.97.10.730 - info@becchisline.com  
www.becchisline.com

Slimming, Lymphatic drainage  
***PRESSOTHERAPY***

***Snell System***

Slimming method  
without operator

---

MADE IN ITALY

www.becchisline.com

---

### Body size and weight

Height: 10 cm

Width: 32 cm

Length: 41 cm

Weight: 5 kg

### Size and overalls set

Height: 12 cm

Width: 46cm

Length: 91 cm

Weight: 6,6 kg

---

### Technical features

Power supply: 230V / 50Hz

Electricity consumption: 350 W

Max pressure: adjustable from 30 to 120 mm hg

Electrical absorption: 160 VA

Security class: II / BF

Sternum protection: 2 fuses 1,6 A

Inflating sectors: 18 independent

Compression: 1 to 15 seconds per sector

General time: 5 to 60 minutes

N° of programs: 4

---

## YOU CAN USE IT AS

### IMPROVEMENT OF CIRCULATION

deflation promotes outflow of extracellular impurities towards the vessels by improving the microcirculation lymphatic

### LYMPHATIC DRAINAGE

has a remarkable draining function, improves both blood circulation and the lymphatic one

### SLIMMING

Localized loss of cm

### DRAINING

Remarkable draining function

